

## ***Integrated Health Care***

18530 Spring Creek Drive, Suite A  
Tinley Park, IL 60477

Phone: 708-468-4070  
Fax: 708-468-4073

[www.ihctech.com](http://www.ihctech.com)

## **Fix for Slow reports in Version 20 Medisoft**

For Medisoft Version 20

*Screen prints are on page 2*

### *Problem or Issue:*

Medisoft Reports are running slower than normal after upgrading to Medisoft Version 20

### **Solution:**

Download then unzip the file v20reportsfix.zip

### **On the Medisoft Server Computer Only:**

Copy the files **repField.adt** and **repField.adi** into the folder ReportDefs that is found inside the Medidata folder on the Medisoft server.

**Then, on each computer** that has Medisoft Version 20 installed on it:

Copy the new file named **FRInit.dll**

Then Paste the file into the “**BIN**” Folder (Replace the existing file).

The file **FRInit.dll** must be copied to every computer that has Medisoft installed.

*The “BIN” Folder is usually located at:*

For Windows 7 Computers: C:\Medisoft\Bin

For Windows XP Computers: C:\Program Files\Medisoft\Bin

You can right Click on the Medisoft Shortcut and go to Properties to verify the Install Location.


Screen prints are on page 2

---

*March 2015*

## Slow reports in Version 20 Medisoft Fix

For Medisoft Version 20



### On the Server Computer Only:

Copy the files  
**repField.adt** and  
**repField.adi**

Copy into the folder **ReportDefs**

This is in the Medidata folder on  
the Medisoft server.

