



**To run the DX Ranking Report:  
Primary Diagnoses Ranking by Code (Charge count)**



In Medisoft click on **Reports - Medisoft Reports**



You can also search for the report – type “diag: then click Find Now

When the new window opens, go to:  
**Plus Pack**  
**Primary Diagnoses Ranking by Code (Charge count)**  
Make sure you use “Charge Count”

**To run the DX Ranking Report:  
Primary Diagnoses Ranking by Code (Charge count)**

Search:  Find Now

Report Name	Folder Name	Date Modified
12. Missing diagnosis codes	Data Analysis	11/11/2009
Diagnoses Ranking (Case count)	Plus Pack	3/8/2013
Diagnoses Ranking (Charge count)	Plus Pack	3/8/2013
Export Diagnosis Data	Focus On Data	
Patient Labels (Last Visit Date & Diagnoses)	Plus Pack	
Patient List by Primary Diagnosis	Plus Pack	
Primary Diagnoses Ranking by Code (Case count)	Plus Pack	
Primary Diagnoses Ranking by Code (Charge count)	Plus Pack	
Primary Diagnoses Ranking by ICD-10 (Case count)	Plus Pack	
Primary Diagnoses Ranking by ICD-10 (Charge count)	Plus Pack	
Primary Diagnoses Ranking by ICD-9 (Case count)	Plus Pack	9/25/2013
Primary Diagnoses Ranking by ICD-9 (Charge count)	Plus Pack	9/25/2013

OK  
Cancel

You can also search for the report – type “Diag:” then click Find Now & choose the report. Make sure you use “Charge Count”

Search

Show all data where the **Date From** is between  and

Show all values of the Date From field.

and the **Attending Provider** is between  and

Show all values of the Attending Provider field.

and the **Date Created** is between  and

Show all values of the Date Created field.

OK Cancel

You can put a date range (such as 1/1/2014 to the present)  
The report will preview to the screen.

**To run the DX Ranking Report:  
Primary Diagnoses Ranking by Code (Charge count)**

**Diagnoses Ranking by Code**

(count based on case primary diagnosis)

Happy Valley Medical Clinic

<b>Rank</b>	<b>Code</b>	<b>Description</b>	<b>Count</b>
<b>1.</b>	<b>847.2</b>	<b>Back Spasm</b>	<b>3</b>
	I9: 847.2	Back Spasm	
	I10:		
<b>2.</b>	<b>110.1</b>	<b>Infected Nail</b>	<b>1</b>
	I9: 110.1	Infected Nail	
	I10:		
<b>3.</b>	<b>110.4</b>	<b>Athlete's Foot</b>	<b>1</b>
	I9: 110.4	Athlete's Foot	
	I10:		
<b>4.</b>	<b>401.9</b>	<b>Hypertension</b>	<b>1</b>
	I9: 401.9	Hypertension	
	I10:		
<b>5.</b>	<b>464.0</b>	<b>Laryngitis, Acute</b>	<b>1</b>
	I9: 464.0	Laryngitis, Acute	
	I10:		
<b>6.</b>	<b>465.9</b>	<b>Upper Respiratory Infection, AC</b>	<b>1</b>
	I9: 465.9	Upper Respiratory Infection, AC	
	I10:		
<b>7.</b>	<b>486</b>	<b>Pneumonia</b>	<b>1</b>
	I9: 486	Pneumonia	
	I10:		
<b>8.</b>	<b>487</b>	<b>Influenza</b>	<b>1</b>
	I9: 487	Influenza	
	I10:		

**Integrated Health Care**

18530 Spring Creek Drive, Suite A, Tinley Park, IL 60477

Phone: 708-468-4070 - Fax: 708-468-4073

Website: [www.ihctech.com](http://www.ihctech.com)